




** PRICES ARE FOB NASHVILLE, TN, FREE SHIPPING ON ALL ORDERS OVER \$1750 **

Picture	Description	Brand - Code	Quantity	UOM	Price	Total
	Pipette 500ml. Autoclavable at 121°C (250°F) for sterilization and decontamination. Disassembly of the instrument is not required, and recalibration is not necessary after autoclaving. Pipettors feature a shorter length for easier handling, especially in hoods. Resistant to UV light.	Eppendorf - 22 47145-7	6	Each	204.98	1,229.89
	Pipetter Multichannel -8/12 Excellent channel-to-channel accuracy and precision. For decontamination or use in sterile applications, detach manifold to autoclave up to 250°F (121°C) without risk to pipetter calibration. Individually replaceable nose cones and user-replaceable seals for easy in-lab maintenance. Easy Calibration* feature allows calibration adjustments in seconds without tools	Brickmann - 050-08-130-3	1	Each	587.00	587.00
	Transferpipettor Ideal for oils, glycerins, solvents, protein solutions, or ether. Handle difficult-to-pipet media such as viscous, foaming, or high vapor pressure liquids with precision and accuracy. Four color-coded models cover range from 100µL to 10mL.	Brickmann - 050-06-100-1	2	Each	264.08	528.16
	Digital Burette II, 50 ML. Buret, Bottletop; Digital; Model 50; Complete with safety valve, telescopic aspiration tube, adjustable discharge tube, three adapters for 40/38/32 mm bottle threads; Capacity: 50 mL	Brickmann - 050-07-510-9	2	Each	737.94	1,475.88
	Volume 1-10ml. Maxitips are designed for use with the Eppendorf Maxipettor repeating pipettor. Maxitip "P" (positive displacement) comes complete with piston and barrel. Recommended for most liquids including viscous solutions. Sterile tips are individually wrapped in RNase-, ATP-, DNase-, pyrogen-, and NACP-free certified Eppendorf Biopur packaging; nonsterile tips are bagged.	Eppendorf - 22291050	6	Case of 100	131.47	788.80
	Safe-Lock Microcentrifuge Tubes, Polypropylene, 1.5 ml, Color red, Ideal for high-speed centrifugation, high-temperature incubations, enzyme assays, cloning procedures. One-piece, polypropylene Safe-Lock tubes have flat-top writing surfaces for easy sample identification. Sealing cap lock snaps closed and stays closed while centrifuging, yet opens easily with one hand. Autoclavable at 121°C (250°F).	Eppendorf - 22363280	5	Case of 500	38.50	192.50
	Microcentrifuge Tubes Description •With attached caps •Capacity 400 µL and 1.5 mL •Disposable polypropylene •Autoclavable at 121 °C •Conical Flex-Tubes™ are available in six colors	Eppendorf - 22363905	4	Case of 1000	56.42	225.68

Picture	Description	Brand - Code	Quantity	UOM	Price	Total
	Combitips 0.5 ML. Combitips pipet tips have polyethylene plungers and polypropylene cylinders. Tips can be used with Repeater, and Repeater Pro models. Nonsterile tips are packaged in bags of 100	Eppendorf - 22495054 (similar to 22266101)	7	Case of 100	74.97	524.79
	Volume 2-200 epTIPS (Eppendorf Totally Integrated Pipetting System) pipet tips satisfy the most stringent requirements of dependability and reproducibility, and perfectly match every Eppendorf pipettor. Tips are also universally applicable for use with pipettors from all other manufacturers. Flexible crown design results in an improved seal with the pipettor and reduced attachment and ejection forces. Tips feature fine volume graduations and a soft, flexible tip cone.	Eppendorf - 22351371 (actual code 022491938)	3	Case of 960	47.36	142.07
	1000 ML filtertips (singletip) Volume 5 -1000. Filter tips are PCR Clean, sterile, have an integrated hydrophobic polyethylene filter, and do not contain "self-sealing" additives. Sample can be easily recovered, even if filter becomes wet as a result of incorrect pipetting.	Eppendorf - 22490346 (similar to 22491253)	5	Case of 96	95.68	478.38
	1000 ML Quality Tips. Volume 5 -1000. Quality tips are precisely molded and feature defined flow dynamics, minimal water retention, high chemical resistance, and high thermal stability. Sterile tips are available in racks or individually wrapped to ensure sterility	Eppendorf - 22350978 (similar to 22491954)	2	Case of 96	51.60	103.21
	Positive Displacement Tips The tip and piston assembly is designed for easy, reliable and secure connection to the pipettor. These disposable positive displacement tips eliminate aerosols and reduce contamination for cleaner, more reproducible sample preparation. Accuracy: ±6 to 2.5%. Precision: <=4.0 to 0.8%. Tips have a 1 to 20µL volume, and are preassembled with pistons and packaged in autoclavable racks.	Eppendorf - 22254159	2	Case of 480 (5 racks of 96)	173.60	347.20

PRICES & AVAILABILITY SUBJECT TO CHANGE WITHOUT NOTICE - WE RESERVE THE RIGHT TO LIMIT QUANTITIES.

*** FREE SHIPPING & PROMOTIONAL OFFERS APPLY ONLY TO "IN STOCK" ITEMS AND ORDERS PAID WITHIN ESTABLISHED TERMS.